

How to use our containers


IT IS STRICTLY PROHIBITED TO PUT ANY OTHER THAN THE SPECIFIED WASTE INTO THE CONTAINERS BELOW. EACH CONTAINER MUST EXCLUSIVELY BE USED FOR THE DISPOSAL OF THE WASTE LISTED ON ITS LABEL.



Used chairside traps/solid collectors, amalgam sludge



Extracted teeth with amalgam fillings



Empty/used amalgam capsules



General safety advice:

- ✓ Always wear protective gloves, a face mask and protective clothing when handling infectious or amalgam waste
- ✓ Never put containers with amalgam waste close to heat sources or in direct sunlight, always keep them as cool as possible
- ✓ Make sure all containers/bags are safely closed before putting them in their slot in the box.

Separator replacement cycle

Separator Cartridge (serial number)	Date of replacement

If you have any questions, just contact us!
Phone: 888-640-2527
E-mail: info@medentex.com
Web: www.medentex.com

Please send all containers back, even if they are not full.

The Silver Service process in short

